

COVID-19 Ag Rapid Saliva Test Device

Cat#: CO-07

For professional in vitro diagnostic use only.

Intended Use

COVID-19 Ag Rapid Saliva Test Device is a rapid chromatographic immunoassay for the qualitative detection of N antigen from SARS-CoV-2 present in human saliva within the first 7 days of symptom onset. This test is for professional used only, as an aid to early diagnosis of SARS-CoV-2 infection in patient.

The result of this test should not be the sole basis for the diagnosis; confirmatory testing is required.

Summary

Coronavirus is a single-stranded positive-sense RNA virus with an envelope of about 80 to 120 nm in diameter. Its genetic material is the largest of all RNA viruses and is an important pathogen of many domestic animals, pets, and human diseases. It can cause a variety of acute and chronic diseases. Common signs of a person infected with a coronavirus include respiratory symptoms, fever, cough, shortness of breath, and dyspnea. In more severe cases, infection can cause pneumonia, severe acute respiratory syndrome, kidney failure, and even death. The 2019 new coronavirus, or "SARS-CoV-2 (COVID-19)" named by the World Health Organization can cause pneumonia epidemic.

The detection results of this kit are for clinical reference only. The result of this test should not be the sole basis for the diagnosis; confirmatory testing is required.

Principle

The COVID-19 Ag Rapid Saliva Test Device uses double antibody sandwich immunoassay. The NC membrane pre-immobilized with monoclonal antibodies against SARS-CoV-2 antigen and anti-mouse polyclonal antibodies, and the colloidal-gold conjugated with monoclonal antibodies specific to SARS-CoV-2 antigen.

If SARS-CoV-2 antigen present in the sample, a complex formed between the anti-SARS-CoV-2 conjugate and the antigen will be caught by the specific anti- SARS-CoV-2 monoclonal coated on the T region. Results appear in 10 to 20 minutes in the form of a red line that develops on the strip.

Whether the sample contains the SARS-CoV-2 antigen or not, the solution continues to migrate to encounter another reagent (an anti-mouse IgG antibody) that binds the remaining conjugates, thereby producing a red line on the region C.

Kit Content

1) Test device (individually packed in a foil pouch with desiccant).

2) Instruction for use.

Materials required but not provided

Timer

Precautions

- For *in vitro* diagnostic use only.
- Do not re-use the test device.
- Do not use after the expiration date.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Perform test at room temperature 15 to 30 °C.

•Wear gloves when hanging the samples, avoid touching the reagent membrane and sample window.

• All samples and used accessories should be treated as infectious and discarded according to local regulations.

• Avoid using blood samples.

Storage and Stability

Store the COVID-19 Ag Rapid Saliva Test Device at 2-30 °C. Do not freeze. All reagents are stable until the expiration dates marked on their outer packaging and buffer vial.

Test Procedure

The COVID-19 Ag Rapid Saliva Test Device is designed for one-step test for COVID-19 antigen in human saliva..

Allow the test device to reach room temperature 15 - 30°C, and instruct the user not to eat, drink, smoke or chew tobacco products for at least 10 minutes prior to use the test.

1) Remove the test device from the sealed pouch and use the device as soon as possible. Best results will be obtained if the assay is performed immediately after opening the foil pouch.

2) Cough deeply twice before collecting the samples.

3) Pull the blue cap off gently by holding the sides to expose the collection pad.

4) Hold the top portion of the device and place the collection pad into the mouth.

5) Rub the collection pad against the cheek and tongue gently in a circular motion about 10 times. And then place the collection pad in the mouth for about $1\sim2$ minutes until the C line show up in the C region.

6) Remove the device from mouth as soon as the C line appears at the C rigion.

7) Place the cap onto the device, lay it on a flat surface.

8) Read results at 10-15minutes after removing device from mouth. Do not read results after 20 minutes.





Fig. A Gently rub the collection pad against each cheek several times.






Fig B. Gently rub the collection pad on top of the tongue.

Fig C. Place the collection pad underneath the tongue.

Interpretation of Results

NEGATIVE:

Only one red band appears in the control region (C), and no band in the test region (T). The negative result indicates that there are no Novel coronavirus antigen in the sample or the number of viral particles is below the detectable range.


POSITIVE:

Two red bands appear. One red band appears in the control region (C), and one red band in the test region (T).

The shade of color may vary, but it should be considered positive whenever there is even a faint band.

INVALID:

No red band appears in the control region (C). The test is invalid even if there is a band on test region (T). Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the test procedure and repeat the test using a new test device.



Do not interpret the result after 20 minutes.



• The COVID-19 Ag Rapid Saliva Test Device is an initial screening test for qualitative detection. Sample collected may contain antigen titles below the reagent's sensitivity threshold, so a negative test result does not exclude infection with novel coronavirus.

• The COVID-19 Ag Rapid Saliva Test Device detects viable and non-viable novel coronavirus antigen. Test performance depends on antigen load in the sample and may not correlate with cell culture performed on the same sample. A positive test does not rule out the possibility that other pathogens may be present. Therefore, the results must be compared with all other available clinical and laboratory information to make an accurate diagnosis.

• A negative test result may occur if the level of extracted antigen in a specimen is below the sensitivity of the test or if a poor quality specimen is obtained.

• Performance of the test has not been established for monitoring antiviral treatment of novel coronavirus.

• Positive test results do not rule out co-infections with other pathogens.

• Negative test results are not intended to rule in other coronavirus infection except the SARS-CoV-2.

• Children tend to shed virus for longer periods of time than adults, which may result in differences in sensitivity between adults and children list.

• A negative result may occur if the concentration of antigen or antibody in a specimen is below the detection limit of the test or if the specimen was collected or transported improperly, therefore a negative test result does not eliminate the possibility of SARS-Cov-2 infection, and should be confirmed by viral culture or a molecular assay or ELISA.

Performance Characteristics

Clinical Evaluation

Clinical evaluation was performed to compare the results obtained by COVID -19 Ag Rapid Saliva Test Device and PCR. The results were summarized below:

Table 1: COVID Ag Rap	d Saliva Test Device vs.	PCR
-----------------------	--------------------------	-----

		COVID Ag Rapid Saliva Test Device		Total Result
		+	-	Kesuit
PCR	+	129	6	135
PCK	-	0	200	200
Total Res	sults	129	206	335

Relative sensitivity: 129/135= 95.6% (95%CI 93.79%~98.66%) Relative specificity: 200/200 >99% (95%CI 97.12%~100%) Overall agreement: (129+200)/(129+0+6+200)×100%=98.21% (95%CI 96.55%~99.36%)

CI: Confidence Interval

Cross Reaction:

The COVID -19 Ag Rapid Saliva Test Device was evaluated with a total of 47 other viruses and bacteria. The results show that the

COVID-19 Ag Rapid Saliva Test Device has no cross-reactivity with other viruses or microorganisms.

Table 2: Cross-reactivity results

Virus/Bacteria/Parasite	Strain	Concentration	Results
MERS-coronavirus	N/A	36 ug/mL	No Cross-Reactivity
	Type 1	1.5E+05TCID50/mL	No Cross-Reactivity
	Type 3	7.5E+05TCID50/mL	No Cross-Reactivity
	Type 5	4.5E+05TCID50/mL	No Cross-Reactivity
	Type 7	1.0E+05TCID50/mL	No Cross-Reactivity
Adenovirus	Type 8	1.0E+05TCID50/mL	No Cross-Reactivity
	Type 11	2.5E+05TCID50/mL	No Cross-Reactivity
	Type 18	2.5E+05TCID50/mL	No Cross-Reactivity
	Type 23	6.0E+05TCID50/mL	No Cross-Reactivity
	Type 55	1.5E+05TCID50/mL	No Cross-Reactivity
	H1N1 Denver	3.0E+07TCID50/mL	No Cross-Reactivity
	H1N1 WS/33	2.0E+07TCID50/mL	No Cross-Reactivity
Influenza A	H1N1 A/Mal/302/54	1.5E+07TCID50/mL	No Cross-Reactivity
	H1N1 New Caledonia	7.6E+07TCID50/mL	No Cross-Reactivity
	H3N2A/Hong Kong/8/68	4.6E+07TCID50/mL	No Cross-Reactivity
	Nevada/03/2011	1.5E+07TCID50/mL	No Cross-Reactivity
Influenza B	B/Lee/40	8.5E+07TCID50/mL	No Cross-Reactivity
	B/Taiwan/2/62	4.0E+07TCID50/mL	No Cross-Reactivity
Respiratory syncytial virus	N/A	2.5E+05TCID50/mL	No Cross-Reactivity
	Bloomington-2	1 × 10 ⁵ PFU/mL	No Cross-Reactivity
Legionella	Los Angeles-1	1 × 10 ⁵ PFU/mL	No Cross-Reactivity
pneumophila	82A3105	$1 \times 10^5 \text{ PFU/mL}$	No Cross-Reactivity
	К	$1 \times 10^5 \text{ PFU/mL}$	No Cross-Reactivity
	Erdman	$1 \times 10^5 \text{ PFU/mL}$	No Cross-Reactivity
Mycobacterium	HN878	$1 \times 10^5 \text{ PFU/mL}$	No Cross-Reactivity
tuberculosis	CDC1551	$1 \times 10^5 \text{ PFU/mL}$	No Cross-Reactivity
	H37Rv	$1 \times 10^5 \text{ PFU/mL}$	No Cross-Reactivity
Streptococcus	4752-98 [Maryland	$1 \times 10^5 \text{ PFU/mL}$	No Cross-Reactivity

(D1)6B-17]		
178 [Poland		No Cross-Reactivity
23F-16]	1 × 10° PFU/mL	
262 [CIP 104340]	$1 \times 10^5 \text{ PFU/mL}$	No Cross-Reactivity
Slovakia 14-10		
[29055]	$1 \times 10^{\circ} \text{ PFU/mL}$	No Cross-Reactivity
Typing strain		
T1[NCIB 11841, SF	$1 \times 10^5 \text{ PFU/mL}$	No Cross-Reactivity
130]		
Mutant 22	1 × 10 ⁵ PFU/mL	No Cross-Reactivity
FH strain of E aton		
Agent	$1 \times 10^5 \text{ PFU/mL}$	No Cross-Reactivity
[NCTC10119]		
36M129-B7	$1 \times 10^5 \text{ PFU/mL}$	No Cross-Reactivity
229E	1.5E+05TCID50/mL	No Cross-Reactivity
OC43	1.5E+05TCID50/mL	No Cross-Reactivity
NL63	1.5E+05TCID50/mL	No Cross-Reactivity
HKU1	1.5E+05TCID50/mL	No Cross-Reactivity
D 0 0000		
Peru2-2002	1.5E+05TCID50/mL	No Cross-Reactivity
IA10-2003	1.5E+05TCID50/mL	No Cross-Reactivity
Type 1	1.5E+05TCID50/mL	No Cross-Reactivity
Type 2	1.5E+05TCID50/mL	No Cross-Reactivity
Type 3	1.5E+05TCID50/mL	No Cross-Reactivity
Type 4A	1.5E+05TCID50/mL	No Cross-Reactivity
N/A	1.5E+05TCID50/mL	No Cross-Reactivity
C-TAN-nCOV		No Corres Desertionit
wuhan strain 01	1.5E+051C1D50/mL	No Cross-Reactivity
	178 [Poland 23F-16] 262 [CIP 104340] Slovakia 14-10 [29055] Typing strain T1[NCIB 11841, SF 130] Mutant 22 FH strain of E aton Agent [NCTC10119] 36M129-B7 229E OC43 NL63 HKU1 Peru2-2002 IA10-2003 Type 1 Type 1 Type 2 Type 3 Type 4A N/A C-TAN-nCOV	178 [Poland $23F-16]$ 1×10^5 PFU/mL 262 [CIP 104340] 1×10^5 PFU/mL Slovakia 14-10 (29055] 1×10^5 PFU/mL Typing strain 1×10^5 PFU/mL TINCIB 11841, SF 1×10^5 PFU/mL 130] 1×10^5 PFU/mL Mutant 22 1×10^5 PFU/mL FH strain of E aton Agent 1×10^5 PFU/mL [NCTC10119] 1×10^5 PFU/mL 229E $1.5E+05TCID50/mL$ OC43 $1.5E+05TCID50/mL$ NL63 $1.5E+05TCID50/mL$ Peru2-2002 $1.5E+05TCID50/mL$ IA10-2003 $1.5E+05TCID50/mL$ Type 1 $1.5E+05TCID50/mL$ Type 1 $1.5E+05TCID50/mL$ Type 1 $1.5E+05TCID50/mL$ Type 3 $1.5E+05TCID50/mL$ Type 4A $1.5E+05TCID50/mL$ N/A $1.5E+05TCID50/mL$

Microbial Interference:

Microbial interference study was performed to evaluate microbial interference effect, using samples spiked at $3 \times \text{LoD}$ SARS-CoV-2 concentration and a high interferent level.

Table 3: Microbial interference Results

No. Microorganism Concentration Results



1	Streptococcus hemolyticus	1×10^5 cfu/ml	No Interference
2	Pseudomonas aeruginosa	1×10^5 cfu/ml	No Interference
3	Staphylococcus aureus	1×10^5 cfu/ml	No Interference
4	Escherichia coli	1×10^5 cfu/ml	No Interference
5	Candida albicans	1×10 ⁵ cfu/ml	No Interference
6	Aspergillus niger	1×10^5 cfu/ml	No Interference

The results show that microorganism listed above has no microbial interference on the negative and positive test results, and these substances do not cross-react with COVID-19 Antigen Rapid Saliva Test Device.

Endogenous Interference:

The COVID -19 Ag Rapid Saliva Test Device was evaluated with a total of 13 endogenous interference substances.

Table 4: Endogenous Interference

Substance	Concentration	Results
Whole Blood	4%	No Interference
Mucin	0.5%	No Interference
Benzocaine	1.5 mg/mL	No Interference
NeilMed	5% v/v	No Interference
CVS Nasal Drops	150//.	No Interference
(Phenylephrine)	15% v/v	No Interference
Oxymetazoline)	15% v/v	No Interference
CVS Nasal Spray	150/ /	
(Cromolyn)	15% v/v	No Interference
Zicam	5% v/v	No Interference
Sore Throat Phenol Spray	15% v/v	No Interference
Tobramycin	4 μg/mL	No Interference
Mupirocin	10 mg/mL	No Interference
Fluticasone Propionate	5% v/v	No Interference
Tamiflu	5 mg/mL	No Interference

The results show that endogenous interference substances listed in above table has no inference effect on the negative and positive test results, and these substances do not cross-react with COVID-19 Antigen Rapid Saliva Test Device.

Food/beverage Interference:

Food/beverage interference study was performed to evaluate the potential interference of food/beverage in saliva samples on the COVID-19 Antigen Rapid Saliva Test Device.

Table 5: Food/beverage interference Result	ts
--	----

Substance	Concentration	Results
Mouth Wash	1%	No Interference
Orange Juice	1%	No Interference
Alcohol	1%	No Interference
MSG	1%	No Interference
Salt	1%	No Interference
Tooth Paste	1%	No Interference
Gum	1%	No Interference
Cough Syrup	1%	No Interference
Sugar	1%	No Interference
Tea	1%	No Interference
Food Color: red	1%	No Interference
Food Color: blue	1%	No Interference
Food Color: green	1%	No Interference
Cranberry Juice	1%	No Interference
Carbonated Cola	1%	No Interference
Baking Soda	1%	No Interference
Cigarette	1%	No Interference

The results show that 1% substance listed in above Table has no inference effect on the negative and positive test results, and these substances have no interference on COVID-19 Antigen Rapid Saliva Test Device.

Manufactured By:

Jiangsu Well Biotech Co., Ltd. No. 9 Changyang Rd., Changzhou, Jiangsu, China 213149 Tel: +86-519-80191660 Fax: +86-519-80191663

European Representative:

Lotus NL B.V. Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands. Email: peter@lotusnl.com Symbols

Symbol	Meaning
Ĩ	Consult instruction for use
IVD	In-Vitro Diagnostic Medical Device
	Manufacturer
LOT	Batch code
\triangle	Caution, consult accompanying documents
*	Keep away from sunlight
\otimes	Do not reuse
X	Temperature Limitation
	Use by date
\sim	Production Date
Σ	Contains sufficient for <n>test</n>
EC REP	Authorized representative in the European Community
CE	Meet the requirements of EC Directive 98/79/EC

Effective Date: 2021.04.11